


# STEP-BY-STEP GUIDE TO USING YOUR SKYRIZI ON-BODY INJECTOR

Skyrizi<sup>®</sup> COMPLETE



Please see **Uses and Important Safety Information** on page 22.

Please see full **Prescribing Information**, including **Medication Guide**, at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>.

  
**Skyrizi**<sup>®</sup>  
risankizumab-rzaa



## USES<sup>1</sup>

SKYRIZI is a prescription medicine used to treat adults with:

- moderate to severe Crohn's disease.
- moderate to severe ulcerative colitis.

## SAFETY CONSIDERATIONS<sup>1</sup>

**Do not use if you are allergic to SKYRIZI. SKYRIZI may cause serious side effects,** including serious allergic reactions, an increased risk of infections, and liver problems.

- Stop using SKYRIZI and get emergency medical help right away if experiencing symptoms of a serious allergic reaction.
- Before starting treatment, your doctor should check for infections and tuberculosis. Tell your doctor right away about any infections or symptoms of infection and about planned or recent vaccines.
- Liver problems may happen while being treated for Crohn's disease or ulcerative colitis that can lead to hospitalization. Your doctor will do liver blood tests before and during treatment and may stop SKYRIZI if liver problems develop.



Please see additional [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf), including [Medication Guide](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>


**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa

# TABLE OF CONTENTS

Use the links below to find the specific information you’re looking for.

- Getting Ready to Inject at Home**..... 4
- Your SKYRIZI On-Body Injector (OBI)**..... 5
- Let the Lights and Sounds Guide You**..... 6
- Step 1.** Get Ready..... 7
- Step 2.** Set Up the OBI..... 11
- Step 3.** Prepare to Inject..... 14
- Step 4.** Inject Medicine ..... 17
- Step 5.** Confirm and Dispose ..... 19
- Remember the 5 P’s**..... 20
- The Skyrizi Complete App** ..... 21



### SEE HOW TO USE THE OBI

Watch Dave give step-by-step instructions at **SKYRIZIOBtraining.com**

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie’s privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa

# GETTING READY TO INJECT AT HOME

Once you've completed your SKYRIZI infusions, you will begin using the SKYRIZI On-Body Injector (OBI)—a device that injects your medicine just under your skin. Your healthcare provider should show you how to properly inject before you start injecting on your own.

Use the **Instructions for Use** that comes with your OBI as your main resource. It will guide you through each step, showing you exactly what to do and how to do it. Then use this guide as a supplemental resource.



## GET ONGOING SUPPORT WITH INJECTING

Your Skyrizi Complete Nurse Ambassador\* can provide support to help you feel prepared for your prescribed SKYRIZI treatment. They can walk you through the OBI injection process, either in-person or virtually, and show you how to inject the right way.

Call **1.866.SKYRIZI** (1.866.759.7494).<sup>†</sup>

SKYRIZI is available in a 180 mg/1.2 mL or 360 mg/2.4 mL subcutaneous injection using the OBI. Your doctor will decide which maintenance dose is best for you.

\*Nurse Ambassadors are provided by AbbVie and do not work under the direction of your health care professional (HCP) or give medical advice. They are trained to direct patients to their HCP for treatment-related advice, including further referrals.

<sup>†</sup>Help is available Monday through Friday, from 8:00 AM to 8:00 PM ET, except for holidays.

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa



# YOUR SKYRIZI OBI


## FRONT

**GRAY DOOR**

**Do not** close gray door without prefilled cartridge inside

**MEDICINE WINDOW**

**PULL TABS**



**ADHESIVE MATERIAL**

**STATUS LIGHT**

**GRAY START BUTTON**  
**Do not** touch until ready to inject

## BACK


**ADHESIVE BACKING**

**NEEDLE COVER**

**NEEDLE INSIDE**  
(under needle cover)

**Do not touch the needle cover or needle**

**CLEAR PLASTIC STRIP**



**LARGE SECTION GREEN PULL TAB**

**SMALL SECTION GREEN PULL TAB**

## SIDE

**NEEDLE COVER**

**NEEDLE INSIDE**  
(under needle cover)

**Do not touch the needle cover or needle**



**DOOR LATCH**  
Opening side has ridges

**Gray door should be slightly open**

**Do not** close gray door without prefilled cartridge inside.

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.


To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE


**Skyrizi**®  
risankizumab-rzaa

# LET THE LIGHTS AND SOUNDS GUIDE YOU


Once you're ready to start your injection, the status light and different sounds of your OBI will guide you through the process. Get to know what each light means to help you stay on track when injecting. Read the **Instructions for Use** before you start injections.




Before you push the large gray start button, your device will have a flashing **blue** light. This means the OBI is On/Ready.




The light will flash **green** after you press the start button. This means the injection has started.



A solid **green** light means the injection has finished.



If the status light shows no light, your OBI is off or finished.



If the status light flashes **red**, it means your OBI is not working properly. Do not continue to use it and call **1.866.SKYRIZI** (1.866.759.7494).



## KEY SYMBOLS USED IN THIS GUIDE



**Pay close attention**



**Caution - Be aware or follow closely**

Follow the tips next to these symbols to help make sure your OBI operates correctly.

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa



# STEP 1: GET READY

Before using the SKYRIZI OBI, you should receive training on how to prepare and inject SKYRIZI from your doctor or a healthcare provider.

Full directions can be found in the **Instructions for Use** that comes with your OBI.

1A



## Check expiration date

- Remove the carton from the refrigerator
- Check expiration date (EXP:) on the carton
- **Do not** use SKYRIZI if the expiration date (EXP:) has passed

1B



- Leave unopened carton at room temperature to warm.
- Keep out of direct sunlight
- Wait at least 45 minutes up to a maximum of 90 minutes to allow SKYRIZI to warm.
  - The OBI will not work if left at **room temperature** for less than **45 minutes**.
- **Do not** warm SKYRIZI in any other way (for example, **do not** warm it in a microwave or in hot water)



## USE THE TIMER IN THE SKYRIZI COMPLETE APP

Make sure SKYRIZI warms to room temperature before injecting. Set a timer for 45-90 minutes in the App.

## Download the App



Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie’s privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa

STEP 1: GET READY (con't.)

1C



Open the carton and remove the plastic tray

- Lift up the flap on the side of the carton
- Take out the plastic tray

1D

Gather all supplies and wash hands

- Place the below items on a clean surface. Wash and dry your hands



Tray with OBI and cartridge



2 alcohol wipes (not included)



Cotton ball or gauze pad (not included)



Sharps disposal container (not included)



Instructions for Use

**NOTE:** Cotton balls, gauze pads, and alcohol wipes are available at your local pharmacy.

Please see [Uses and Important Safety Information](#) on page 22.


Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie’s privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>



## STEP 1: GET READY (con't.)

1E



### Remove the white paper tray seal

- Locate the black arrow
- Peel away the white paper tray seal from the plastic tray
- **Do not** use the OBI and prefilled cartridge if the white paper tray seal is missing or damaged. Call the pharmacy where you obtained your OBI for return instructions

1F



### Lift the plastic cover

- Locate the rounded opening on the top cover
- Insert your pointer finger in the opening and place your thumb on the opposite side
- Lift the cover to remove and set aside

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.


To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa

## STEP 1: GET READY (con't.)

1G



### Inspect the OBI

- Check that the OBI is intact and undamaged
- **Do not** close the gray door before the prefilled cartridge is loaded
- Hold the OBI with the gray door at the top. The gray door should be slightly open
- If the gray door does not open, press in firmly on the gray door ridges (left side of the door) and swing open the door
- **Do not** use OBI if you drop it, discover missing pieces, or it is damaged
- **Do not** let the OBI get wet with water or any other liquids



START BUTTON



NEEDLE COVER



**Do not** touch the gray start button until it is time to inject. It can be pressed 1-time only.

If the gray start button is pressed before placing it on your body, the OBI can no longer be used. Call **1.866.SKYRIZI** (1.866.759.7494) for help.

**Do not** touch the needle cover or the needle.



Proceed to STEP 2: **SET UP THE OBI**

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa



## STEP 2: SET UP THE OBI

2A




### Fully open the gray door



Avoid touching the needle cover or the needle on the back of the OBI. The needle is behind the needle cover.

- Swing the gray door all the way to the right to open it
- If the gray door does not open, press in firmly on the gray door ridges (left side of the door) and swing open the door
- **Do not** close the gray door before the prefilled cartridge is loaded
- **Do not** touch the needle cover or needle
- Set the OBI aside with the door open while you inspect the prefilled cartridge

2B



### Inspect the prefilled cartridge

- Carefully remove the prefilled cartridge from the plastic tray
- **Do not** twist or remove cartridge top
- Check the prefilled cartridge to make sure:
  - The expiration date (EXP:) has not passed.
  - The liquid should look clear to yellow and may contain tiny white or clear particles. It is normal to see one or more bubbles
  - The prefilled cartridge parts appear intact, and the clear plastic is not cracked or broken
- **Do not** use if the expiration date (EXP:) has passed
- **Do not** use if the liquid is cloudy, discolored, or contains flakes or large particles
- **Do not** use if the liquid has been frozen (even if thawed)
- **Do not** use the prefilled cartridge if you drop it, discover missing pieces, or it is damaged

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa



## STEP 2: SET UP THE OBI (con't.)

2C

**Clean the smaller bottom tip of the prefilled cartridge**

- Locate the smaller bottom tip of the prefilled cartridge
- Clean smaller bottom tip of the prefilled cartridge with an alcohol wipe. Make sure to use the alcohol wipe to clean the center of the smaller bottom tip of the prefilled cartridge
- **Do not** touch the smaller bottom tip of the prefilled cartridge after cleaning

2D

**Insert the cleaned prefilled cartridge straight into the OBI**

- **Do not** twist or remove the prefilled cartridge top
- Insert the smaller bottom tip of the prefilled cartridge into the OBI first



- **Firmly push down on the prefilled cartridge top until you hear a "CLICK"**
- After loading the prefilled cartridge, you may see a few drops of medicine on the back of the OBI. This is normal



**Make sure that you start the injection within 5 minutes after inserting the prefilled cartridge into the OBI. Waiting will dry out the medicine and the OBI will not work afterwards.**



Remember, you can watch step-by-step instructions anytime at **SKYRIZIOBtraining.com**. If you have any questions, reach out to your Nurse Ambassador at **1.866.SKYRIZI** (1.866.759.7494).

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa




# STEP 2: SET UP THE OBI (con't.)

2E



CLOSE THE GRAY DOOR FIRMLY, UNTIL YOU HEAR A “SNAP”



## Close the gray door by swinging it to the left and squeezing firmly until you hear a “SNAP”

- Swing the gray door to the left, then squeeze firmly and listen for the gray door to “snap” shut
- The gray door should stay locked after loading the prefilled cartridge
- **Do not** close the gray door if the prefilled cartridge is not fully inserted or is missing
- Proceed without delay to Prepare to Inject

Set a timer and start the injection within 5 minutes after loading the prefilled cartridge into the OBI.



Proceed to STEP 3: **PREPARE TO INJECT**

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie’s privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa

## STEP 3: **PREPARE TO INJECT**

Keep electronic devices including cell phones at least 12 inches (30 cm) from the OBI until injection is complete

3A

**STOMACH****THIGH**

### Pick your injection area

Locate a firm and flat area on your bare skin. Make sure the area is clear of clothing:

- Your abdomen (belly) at least 2 inches from your navel (belly button).
- The front of your left thigh or your right thigh
- **Do not** inject into areas of the skin that naturally fold or bulge because the OBI could fall off during wear
- **Do not** inject into skin that is sore, bruised, red, hard, scarred, or has stretch marks, moles, or excessive hair. You can trim excess hair from the injection area for easier access, if needed

3B



### Clean the injection area

- Clean the injection area with an alcohol wipe
- Allow the injection area to completely dry
- **Do not** touch the cleaned injection area before placing the OBI on the skin

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa



# STEP 3: PREPARE TO INJECT (con't.)



3C



LARGE SECTION



SMALL SECTION

## Peel both tabs to expose adhesive

- Turn the OBI over to find both **green** pull tabs



Avoid touching the needle cover (needle inside)

- Peel away the large section using the **green** pull tab to expose the adhesive
- Peel away the small section using the **green** pull tab to expose the adhesive. This will remove the clear plastic strip, activating the OBI
- **Do not** pull the adhesive material off the OBI or allow the sticky side to fold and stick to itself



STATUS LIGHT  
FLASHING BLUE



- Check the status light when the OBI beeps
- The status light will flash **blue** when the OBI is activated
- If the status light does not flash **blue**, call **1.866.SKYRIZI** (1.866.759.7494)
- **Do not** press the gray start button yet
- **Do not** touch the needle cover or the needle



If the status light flashes red, and beeps are heard, the OBI is not working properly. Do not continue to use it.

Call **1.866.SKYRIZI** (1.866.759.7494) for assistance. If OBI is attached to the body, carefully remove from the skin.



beep-beep-beep-beep-beep

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie’s privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa

## STEP 3: PREPARE TO INJECT (con't.)

3D



STOMACH




THIGH

### Prepare the OBI for placement

- For the stomach, move and hold the skin to create a firm, flat surface for injection at least 2 inches from your navel (belly button). Make sure to adjust your posture (sit up straight) to avoid skin folds and bulges
- You do not need to pull the skin flat for the front of left thigh or right thigh
- Make sure to position the OBI so that you can see the **blue** status light

3E



### Place the OBI on your skin and press firmly

- When the **blue** light flashes, the OBI is ready. Place the OBI onto the cleaned skin with the status light visible
- **Do not** place the OBI on clothes. Only place on bare skin
- Run your finger around the adhesive material to secure it. Make sure all of the adhesive is attached to your skin
- **Do not** move or adjust the OBI after it has been placed on your skin



When the OBI is on your skin, proceed to  
STEP 4: **INJECT MEDICINE**

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa



## STEP 4: INJECT MEDICINE

4A



FIRMLY PRESS  
START BUTTON  
**ONE TIME ONLY**  
UNTIL YOU HEAR  
A “CLICK.”



LISTEN FOR **BEEP-BEEP-BEEP**  
THEN CHECK FOR FLASHING  
GREEN STATUS LIGHT



### Press the gray start button firmly to start the injection

- **Firmly** press the gray start button until you hear a click. Then release the gray start button



**Do not** touch the gray start button until you place the OBI loaded with the prefilled cartridge onto your skin and are ready to inject. You can only press the start button 1-time. If the gray start button is pressed before the OBI is placed on your body, it can no longer be used.

- You may feel a needle pinch
- Check the status light when the OBI beeps
- After starting the injection, the status light will continuously flash **green**
- You will hear pumping sounds as the OBI delivers the medicine



**Do not** continue to use the OBI if status light flashes **red** and beeps. Carefully remove from skin and call **1.866.SKYRIZI** (1.866.759.7494).

4B



IT MAY TAKE UP  
TO **5 MINUTES**  
TO COMPLETE  
THE FULL DOSE OF  
SKYRIZI MEDICINE



### Wait for the injection to finish

- It may take up to 5 minutes to complete the full dose of SKYRIZI medicine.
- During the injection:
  - The status light will continue to flash **green**
  - You will hear pumping sounds as the OBI continues delivering the medicine
  - Moderate physical activities can be done, such as walking, reaching, and bending
- The OBI will automatically stop when the injection is finished



**Do not** continue to use the OBI if status light flashes **red** and beeps. Carefully remove from skin and call **1.866.SKYRIZI** (1.866.759.7494).

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>


**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa



## STEP 4: INJECT MEDICINE (con't.)

4C



### Injection is finished

- The OBI will stop on its own
- After the injection is finished, you will hear beeps and the status light will change to solid **green**. Check the status light. If it has changed to solid **green**, this means that the injection is finished
- When the injection is finished, proceed to next step

4D



### Remove the OBI

- **Do not** put your fingers on the back side of the OBI when removing it from your skin
- When the injection is done, grab the corner of the adhesive to carefully peel the OBI from the skin



Avoid touching the needle cover or needle on the back of the OBI

- After removing the OBI, you will hear several beeps and the status light will turn off
- The needle cover will cover the needle when the OBI is removed from the skin
- It is normal to see a few small drops of liquid on your skin after removing the OBI. If you see more than a few small drops of liquid left on your skin, call **1.866.SKYRIZI** (1.866.759.7494)
- Press a cotton ball or gauze pad over the injection site on your skin and hold for 10 seconds
- **Do not** rub the injection site
- Slight bleeding at the injection site is normal



Proceed to STEP 5: **CONFIRM AND DISPOSE**

Please see Uses and Important Safety Information on page 22.

Please see full Prescribing Information, including Medication Guide, at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa



## STEP 5: CONFIRM AND DISPOSE

5A



### Check the OBI

- Inspect the medicine window and status light. Check to see that the used white plunger fills the entire medicine window and the solid **green** light turns off, letting you know that all the medicine has been injected



If the used white plunger does not fill the window, call **1.866.SKYRIZI** (1.866.759.7494).

5B



### Put the used OBI in an FDA-cleared sharps disposal container right away after use

- Leave the prefilled cartridge in the OBI
- **Do not** throw away (dispose of) the used OBI in your household trash
- For more information, see the “**Disposing of SKYRIZI**” section of the **Instructions for Use**



**Great job!** Pat yourself on the back for completing your injection. Remember, you can do this!

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie’s privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**SKYRIZI**® COMPLETE

**SKYRIZI**®  
risankizumab-rzaa



# REMEMBER THE 5 P's

1




## PREPARE

Prepare to inject by gathering alcohol wipes, a gauze pad or cotton ball, and your sharps disposal container, while your SKYRIZI is warming to room temperature outside of the fridge. Wash and dry your hands.



**Make sure that you start the injection within 5 minutes after inserting the prefilled cartridge into the OBI. Waiting will dry out the medicine and the OBI will not work afterwards.**

2



## PICK

Pick the place on your body where you're going to inject.

3



## PEEL

Peel the 2 **green** pull tabs to reveal the adhesive without touching the needle.

4



## PLACE

Place the OBI onto your skin. Make sure you can see the status light.

5



## PRESS

Press and release the gray start button to begin. You'll see a flashing **green** light. When the injection is complete, the **green** light goes solid.



### GET 1-TO-1 SUPPORT WITH INJECTING

Get injection training from your Nurse Ambassador, in person, by phone, or virtually. Call **1.866.SKYRIZI** (1.866.759.7494). Or get answers anytime with **24/7 Live Chat** through the **Skyrizi Complete App**.

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi** COMPLETE

**Skyrizi**  
risankizumab-rzaa






## KEEP TRACK OF INJECTIONS WITH THE SKYRIZI COMPLETE APP

The App puts helpful resources right in the palm of your hand.

- ✓ **Watch** infusion and injection videos to see what to expect
- ✓ **Submit rebates.\*** If eligible, you can submit claims if you have an outstanding infusion bill or have already paid out of pocket for SKYRIZI OBI prescriptions and lab tests
- ✓ **Find** resources that can help you build a routine so you feel prepared switching from infusions to the OBI
- ✓ **Get** answers anytime with 24/7 Live Chat
- ✓ **Access** your Sskyrizi Complete Savings Card, if eligible, and add it to your mobile wallet so you're always ready to pay for treatment



### Download the App



Download on the  
**App Store**



GET IT ON  
**Google Play**

\*For eligible, commercially insured patients only. Please see Terms and Conditions on page 22.

Please see [Uses and Important Safety Information](#) on page 22.

Please see full [Prescribing Information](#), including [Medication Guide](#), at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor.

To learn more about AbbVie's privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>

**Skyrizi**® COMPLETE

**Skyrizi**®  
risankizumab-rzaa

# USES AND IMPORTANT SAFETY INFORMATION

## ABOUT SKYRIZI® (risankizumab-rzaa)<sup>1</sup>

### SKYRIZI USES<sup>1</sup>

SKYRIZI is a prescription medicine used to treat adults with:

- moderate to severe Crohn’s disease.
- moderate to severe ulcerative colitis.

### IMPORTANT SAFETY INFORMATION<sup>1</sup>

#### What is the most important information I should know about SKYRIZI® (risankizumab-rzaa)?

**SKYRIZI is a prescription medicine that may cause serious side effects, including:**

#### Serious allergic reactions:

- Stop using SKYRIZI and get emergency medical help right away if you get any of the following symptoms of a serious allergic reaction:
  - fainting, dizziness, feeling lightheaded (low blood pressure)
  - swelling of your face, eyelids, lips, mouth, tongue, or throat
  - trouble breathing or throat tightness
  - chest tightness
  - skin rash, hives
  - itching

#### Infections:

SKYRIZI may lower the ability of your immune system to fight infections and may increase your risk of infections. Your healthcare provider should check you for infections and tuberculosis (TB) before starting treatment with SKYRIZI and may treat you for TB before you begin treatment with SKYRIZI if you have a history of TB or have active TB. Your healthcare provider should watch you closely for signs and symptoms of TB during and after treatment with SKYRIZI.

- Tell your healthcare provider right away if you have an infection or have symptoms of an infection, including:
  - fever, sweats, or chills
  - cough
  - shortness of breath
  - blood in your mucus (phlegm)
  - muscle aches
  - warm, red, or painful skin or sores on your body different from your psoriasis
  - weight loss
  - diarrhea or stomach pain
  - burning when you urinate or urinating more often than normal

**Do not use SKYRIZI if you are** allergic to risankizumab-rzaa or any of the ingredients in SKYRIZI. See the Medication Guide or Consumer Brief Summary for a complete list of ingredients.

#### Before using SKYRIZI, tell your healthcare provider about all of your medical conditions, including if you:

- have any of the conditions or symptoms listed in the section “What is the most important information I should know about SKYRIZI?”
- have an infection that does not go away or that keeps coming back.
- have TB or have been in close contact with someone with TB.
- have recently received or are scheduled to receive an immunization (vaccine). Medicines that interact with the immune system may increase your risk of getting an infection after receiving live vaccines. You should avoid receiving live vaccines right before, during, or right after treatment with SKYRIZI. Tell your healthcare provider that you are taking SKYRIZI before receiving a vaccine.
- are pregnant or plan to become pregnant. It is not known if SKYRIZI can harm your unborn baby.
- are breastfeeding or plan to breastfeed. It is not known if SKYRIZI passes into your breast milk.
- become pregnant while taking SKYRIZI. You are encouraged to enroll in the Pregnancy Registry, which is used to collect information about the health of you and your baby. Talk to your healthcare provider or call 1-877-302-2161 to enroll in this registry.

**Tell your healthcare provider about all the medicines you take,** including prescription and over-the-counter medicines, vitamins, and herbal supplements.

#### What are the possible side effects of SKYRIZI?

**SKYRIZI may cause serious side effects. See “What is the most important information I should know about SKYRIZI?”**

**Liver problems may happen while being treated for Crohn’s disease or ulcerative colitis:** A person with Crohn’s disease who received SKYRIZI through a vein in the arm developed changes in liver blood tests with a rash that led to hospitalization. Your healthcare provider will do blood tests to check your liver before, during, and at least up to 12 weeks of treatment, and may stop treatment with SKYRIZI if you develop liver problems. Tell your healthcare provider right away if you notice any of the following symptoms: unexplained rash, nausea, vomiting, stomach (abdominal) pain, tiredness (fatigue), loss of appetite, yellowing of the skin and eyes (jaundice), and dark urine.

**The most common side effects of SKYRIZI in people treated for Crohn’s disease and ulcerative colitis include:** upper respiratory infections, headache, joint pain, stomach (abdominal) pain, injection site reactions, low red blood cells (anemia), fever, back pain, urinary tract infection, and rash.

These are not all the possible side effects of SKYRIZI. Call your doctor for medical advice about side effects.

Use SKYRIZI exactly as your healthcare provider tells you to use it. SKYRIZI (risankizumab-rzaa) is available in a 600 mg/10 mL vial for intravenous infusion and a 180 mg/1.2 mL or 360 mg/2.4 mL single-dose prefilled cartridge with on-body injector.

**You are encouraged to report negative side effects of prescription drugs to the FDA. Visit [www.fda.gov/medwatch](http://www.fda.gov/medwatch) or call 1-800-FDA-1088.**

**If you are having difficulty paying for your medicine, AbbVie may be able to help. Visit [AbbVie.com/PatientAccessSupport](http://AbbVie.com/PatientAccessSupport) to learn more.**

**Reference:** 1. SKYRIZI [package insert]. North Chicago, IL: AbbVie Inc.

#### Skyrizi Complete Savings Card Terms and Conditions

Eligibility: Available to patients with commercial insurance coverage for SKYRIZI® (risankizumab-rzaa) who meet eligibility criteria. This co-pay assistance program is not available to patients receiving prescription reimbursement under any federal, state, or government-funded insurance programs (for example, Medicare [including Part D], Medicare Advantage, Medigap, Medicaid, TRICARE, Department of Defense, or Veterans Affairs programs) or where prohibited by law. Offer subject to change or termination without notice. Restrictions, including monthly maximums, may apply. This is not health insurance. **For full Terms and Conditions, visit [www.skyrizi.com/savings-card-terms](http://www.skyrizi.com/savings-card-terms) or call 1.866.SKYRIZI for additional information. To learn about AbbVie’s privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>.**

US-SKZG-240351

**Please see full Prescribing Information, including Medication Guide, at [www.rxabbvie.com/pdf/skyrizi\\_pi.pdf](http://www.rxabbvie.com/pdf/skyrizi_pi.pdf) and discuss with your doctor. To learn more about AbbVie’s privacy practices and your privacy choices, visit <https://abbv.ie/corpprivacy>**



© 2024 AbbVie. All rights reserved. SKYRIZI® and its design are registered trademarks of AbbVie Biotechnology Ltd.